

## soft air® simplex wds

Stage II


Dynamic ulcer pressure relief system with large cells for prevention of ulcers and therapy support for low to medium risk/stage II



ADL® GmbH  
 An der Germania Brauerei 1  
 48159 Münster  
 Germany  
 Telefon: +49 (0) 251 149896-0  
 Telefax: +49 (0) 251 149896-99

ADL® GmbH  
 Produktionsstätte  
 und Servicecenter  
 Im Flürchen 1  
 99334 Amt Wachsenburg  
 Germany  
 Telefon: +49 (0) 36202 7519-0  
 Telefax: +49 (0) 36202 7519-12

E-Mail: [info@adl-gmbh.de](mailto:info@adl-gmbh.de)  
 Internet: [www.adl-gmbh.de](http://www.adl-gmbh.de)



### Technical details

	soft air® simplex wds
Art. No.	110 470-wds
Material	Mattress: polyvinyl · Pump: ABS housing · Cover: PU
Warranty	2 years
Electrical values	220 – 240 V · 50 – 60 Hz/7 W
Fuse	500 mA/slow
Dimensions (inflated)	196 x 88 x 12,5 cm (77" x 34½" x 5")
Weight	Mattress: ca. 5,6 kg (12 lbs) · Pump: ca. 1,3 kg (3 lbs) Cover: ca. 0,9 kg (2 lbs)
Cleaning	Thermal disinfection
Cover	Washable at up to 95 °C, dryable at up to 50 °C
Recommended patient weight	40 – 130 kg (5 – 21 st)
Application area	Ulcer/Decubitus prevention and therapy support from low to medium risk/stage II